

Alternate Design Configurations

A wearer may prefer one progressive lens over another based on their expectations and visual demands. **Ascent Steady Progressive** lenses can be configured for patients' lifestyles.

Ascent Steady Progressive | Near Vision



Ideal wearer: expert or novice progressive lens wearers who have a need for a high value lens and visual demands that call for a larger more usable near area.

Ascent Steady Progressive | Distance Vision



Ideal wearer: expert or novice progressive lens wearers who have a need for a high value lens and visual demands that call for a larger more usable distance area.

Ascent Steady Progressive | Intermediate Vision



Ideal wearer: first-time wearers, wearers who have experienced non-adapt with other progressive lenses, and those who have a need for a high value lens and visual demands that call for a larger more usable intermediate area.

— Ascent Steady Progressive standard configuration



Ascent Steady Progressive Lenses

Free-form progressive lens
with high wearer satisfaction
in an entry-level technology.

Ascent Steady Progressive

A great product portfolio meets all market needs and helps create differentiation. Schroeder Optical offers a non-personalized lens that is entirely different than any other in its market segment. Alternative design configurations are available for specific visual needs or lifestyles.

Ascent Steady Progressive lenses incorporate **Steady Methodology** to enhance image stability and reduce swim effect. Swim effect is a non-stable perception of the surroundings that causes discomfort, and reduces overall lens satisfaction, making adaptation more difficult. **Ascent Steady Progressive** lenses provide wearers with more stable and natural vision.

Superior image stability for more comfortable vision

Steady Methodology is a technological breakthrough in free-form, digital lenses. In addition to controlling for unwanted cylinder power, Steady Methodology addresses unwanted changes to mean power in the lateral areas of the lens. This improves peripheral visual acuity, reduces swim effect, and provides superior image stability and offers more comfortable vision.

Ideal Wearer

Those in need of a high value, reliable solution. Part time eyeglass wearers and wearers with lower prescription and addition powers.



Ascent Steady Progressive | Standard Configuration

Features & Benefits

- Non-compensated progressive lens
- Near, distance, & intermediate enhancements available
- Technology: Steady Methodology
- Precise focus for all working distances
- Higher image stability for reduced swim effect
- Improvement of peripheral visual quality in the distance zone
- High-value solution
- Compatible with virtually any lens material & coating
- MFHs: 14, 15, 16, 17, & 18 mm or variable